

Barcode Scanner — Setup Card

ADEYSONS HEALTH LLC · COMMUNITY HEALTH WORKERS | PHLEBOTOMY | MEDICAL COURIER · CLIA
#45D2344227

Keep this in the vehicle. Works with any Bluetooth barcode scanner.

One-time setup (about 5 minutes)

- 1 Set the scanner to **HID / Keyboard mode** (most ship this way). If not, scan the “Bluetooth HID / keyboard” barcode in the scanner’s manual.
- 2 On your phone open **Settings** → **Bluetooth**. Put the scanner in pairing mode (a button, or a “pair” barcode in the manual).
- 3 Tap the scanner in the phone’s Bluetooth list until it says **Connected**. It now works like a wireless keyboard.

Every shift (10 seconds)

- 1 Open **adeysonshealth.com/scan** in your phone browser (add it to your home screen).
- 2 Enter your **PIN**, choose the **status** (Picked up / In transit / Delivered / Received), and pick or start a **Trip**.
- 3 Tap **Bluetooth scanner**, then tap once inside the “Scan into this box” field.
- 4 **Pull the trigger on each tube, one after another**. Every scan logs instantly — no button to press between tubes.
- 5 At the lab, switch status to **Delivered** (or **Received**) and scan again, then tap **Reconcile trip** — you’ll get an alert if any tube is missing.

Good to know

- 1 A **2D / QR scanner** reads a tube’s QR and advances **that tube** (best for tracking).
- 2 A **1D barcode scanner** reads the order barcode and advances the **whole order** at once.
- 3 The box accepts scans **with or without** an Enter key — any brand works, worldwide.
- 4 No scanner? Use **Camera mode** on the same screen, or your phone’s camera on the label QR.

Recommended hardware: a Bluetooth 2D “ring” or handheld scanner (search “Bluetooth 2D QR ring scanner”, ~\$25–45). Any HID model works.

Generated 09/28/2026 · Questions? info@adeysonshealth.com